


PRODUCTS FOR HYGIENE AND DISINFECTION OF DENTAL SURGERIES



PULI-JET PLUS NEW	
Item Code	060900
Quantity	4 PK
Concentrated disinfectant-detergent for aspirators	
Application: SINGLE-SURGERY AND CENTRALISED ASPIRATION SYSTEMS	
When to use: At the end of the working day, after surgical interventions on high infection risk patients, for a detergent-disinfectant wash.	
How to use:	
<ul style="list-style-type: none"> For a cleaning-disinfectant action: prepare a solution at 0,8% in hot water (about +50°C). For a cleaning-sanitising action: prepare a solution at 0,4% in hot water (about +50°C). Suck from the Pulse Cleaner through the equipment. 	
With 1 litre of product it is possible to obtain 250 litres of sanitising solution at 0,4% and 125 litres of disinfectant solution at 0,8%. Let it work during the night or for at least 1 hour.	

Range of action:
<ul style="list-style-type: none"> Detergent Anti-scale agent (indispensable for wet ring pumps) Disinfectant: – bactericide – fungicide – tuberculicide – Virucidal (HIV, HBV, HCV)

Characteristics: It does not contain phenol. Orange-scented.

Precautions: Use protective gloves when preparing the solution.

Certificates:
<ul style="list-style-type: none"> CE 0434 Class IIa-Medical Device Obtained according to the following applicable rules: – UNI EN 13727-14561 – UNI EN 13624-14562 – Pr EN 14563 – UNI EN 14476

PULI-JET CLASSIC	
Item Code	040715
Quantity	4 PK
Sanitising detergent for the cleaning of aspirators	
Application: SINGLE-SURGERY AND CENTRALISED SUCTION SYSTEMS	
When to use: Midworking day, after surgical interventions, for a detergent-sanitising wash.	
How to use:	
<ul style="list-style-type: none"> Prepare a solution at 5% in hot water (about +50°C). Suck from the Pulse Cleaner through the equipment. 	
With 1 litre of product it is possible to obtain 20 litres of solution at 5%.	
Range of action:	
<ul style="list-style-type: none"> Detergent Sanitising-antibacterial (let it work during the night or for at least 1 hour). 	
Precautions: Use protective gloves when preparing the solution.	

PULSE CLEANER	
Item Code	040720
Quantity	SINGLE
You can obtain a thorough cleaning of the aspiration system with:	
<ul style="list-style-type: none"> “Pulse Cleaner” creating the necessary turbulence. “Puli-Jet Plus” disinfectant-detergent. Antifoam tablets, inhibiting foam formation and reducing the bacterial rate. 	
Pulse Cleaner™ is a new turbulence-maker. To prepare the Puli-Jet Plus solution, pour the concentrated product in warm water (+50°C) and then place the terminals on the pistons.	
The cleansing aspiration will start without any further operation.	
The solution left in the container can be used to sanitise the spittoon bowl.	
Do not rinse the aspirator; do not turn on the spittoon-rinsing tap: Puli-Jet Plus will go on sanitising and dissolving blood and mucus all night long.	

ANTIFOAMING DISINFECTANT TABLETS	
Item Code	040827
Quantity	SINGLE
Disinfectant antifoam tablets for dental aspirators	
Application: SINGLE-SURGERY AND CENTRALISED SUCTION SYSTEMS	
When to use: Every day, after cleaning and disinfection.	
How to use: At the end of every working day, after cleaning, place 1 or 2 tablets in the first filter: when aspirated, the liquid will take up the right quantity of product.	
Range of action: 1 tablet is soluble within 1 average working day.	
<ul style="list-style-type: none"> Antifoam Disinfectant: – bactericide – candidicide 	
Precautions:	
<ul style="list-style-type: none"> Do not remove the protective film: it dissolves in the aspirated liquid. For small filters: remove the film, then break the tablet in half by pressing along the pre-cut notch. Use protective gloves. 	
Certificates:	
<ul style="list-style-type: none"> CE 0434 Class IIa-Medical Device Reference rules: – UNI EN 13727 – UNI EN 13624 	

ECO-JET 1 SPRAY	
Product Code	040745 Size 500ml Qty 4 PK
Product Code	040764 Size 5 litre Qty SINGLE
High-level surface disinfectant-detergent with didecyl ammonium chloride	
Application: IRREGULAR, ROUGH OR DIFFICULT TO REACH SURFACES	
When to use: After every intervention.	
How to use: Ready-to-use.	
<ul style="list-style-type: none"> Spray onto the surface from 20/30 cm. Let it work for 5 minutes. Wipe with a dry sterile cloth. 	
Range of action:	
<ul style="list-style-type: none"> Detergent Disinfectant: bactericide, fungicide, virucide (HBV, HCV, HIV) and tuberculocide (proven effective within less than 15 minutes). 	
Characteristics:	
<ul style="list-style-type: none"> Cool, flower-scented perfume. It contains less than 10% alcohol and is particularly mild on the chair covering material. 	
Precautions:	
<ul style="list-style-type: none"> It is generally safe on hands, it can be used without wearing gloves. It has been tested on many materials. As a precaution, we recommend carrying out a preliminary test on a part which is not in sight. 	
Certificate: 11a Class Medical Device	

Reference standards:
<ul style="list-style-type: none"> CE 0434 UNI EN 14476 (HIV, HBV, HCV) UNI EN 13697 UNI EN 14348 Skin Patch Test

LUBRI-JET SILICONE LUBRICANT DROPS	
Item Code	040727
Quantity	-
Silicone lubricant for the periodic lubrication of the hinge of pincers, scissors and other instruments of the dental surgery	
Application: SURGICAL INSTRUMENTS with parts subject to friction (scissors, pincers, etc)	
When to use: When needed.	
How to use: Dribble a drop on the parts under friction. To be used on clean and dry instruments before sterilisation.	
Range of action: Lubrication	
Characteristics:	
<ul style="list-style-type: none"> It contains high-quality silicone oil. Dust-repellent. 	
Precautions:	
<ul style="list-style-type: none"> No particular storage precaution is required. Avoid contact with eyes. 	

HYGIENE AND DISINFECTION PROGRAM FOR DENTAL ASPIRATORS

PRODUCT	APPLICATION	WHEN TO USE	HOW TO USE	RANGE OF ACTION	CERTIFICATES AND REF. RULES	PRECAUTIONS
PULI-JET PLUS NEW – Concentrated disinfectant-detergent for aspirators						
	Single-surgery and centralised aspiration systems.	At the end of the working day, after surgical interventions on high infection risk patients, for a detergent-disinfectant wash.	For a cleaning-disinfectant action: prepare a solution at 0,8% in hot water (about +50°C). For a cleaning-sanitising action: prepare a solution at 0,4% in hot water (about +50°C). Suck from the Pulse Cleaner through the equipment. With 1 litre of product it is possible to obtain 250 litres of sanitising solution at 0,4% and 125 litres of disinfectant solution at 0,8%. Let it work during the night or for at least 1 hour.	Detergent Anti-scale agent (indispensable for wet ring pumps) Disinfectant: • bactericide • fungicide • tuberculicide • Virucidal (HIV, HBV, HCV)	CE 0434 Class IIa-Medical Device UNI EN 13727-14561 UNI EN 13624-14562 Pr EN 14563 UNI EN 14476	Use protective gloves when preparing the solution.
PULI-JET CLASSIC – Sanitising detergent for aspirators						
	Single-surgery and centralised aspiration systems.	Midworking day, after surgical interventions, for a detergent-sanitising wash.	Prepare a solution at 5% in hot water (about +50°C). Suck from the Pulse Cleaner through the equipment. With 1 litre of Puli-Jet Classic it is possible to obtain 20 litres of solution at 5%. 	Detergent Sanitising-antibacterial		Use protective gloves when preparing the solution.
DISINFECTANT ANTIFOAM TABLETS – Disinfectant antifoam tablets for aspirators						
	Single-surgery and centralised aspiration systems.	At the end of every working day, after cleaning and disinfection.	At the end of every working day, after cleaning, place 1 or 2 tablets in the first filter: when aspirated, the liquid will take up the right quantity of product.	Antifoam Disinfectant: • bactericide, • candidicide	CE 0434 Class IIa-Medical Device UNI EN 13727 UNI EN 13624	Do not remove the protective film: it dissolves in the aspirated liquid. For small filters: remove the film, then break the tablet in half by pressing along the pre-cut notch. Use protective gloves.
ECO-JET 1 SPRAY – High-level surface disinfectant-detergent						
	Irregular, rough or difficult to reach surfaces.	After every intervention.	Ready-to-use. • Spray onto the surface from 20/30cm. • Let it soak for 5 minutes. • Wipe with a dry sterile cloth.	Detergent Disinfectant: bactericide, fungicide, virucide (HBV, HCV, HIV) and tuberculocide (proven effective within less than 15 minutes).	CE 0434 UNI EN 14476 (HBV, HCV, HIV) UNI EN 13697 UNI EN 14348 Skin Patch Test	It is generally safe on hands, it can be used without wearing gloves. It has been tested on many materials. As a precaution, we recommend carrying out a preliminary test on a part which is not in sight.
LUBRI-JET DROP – Lubricant						
	Surfaces exposed to rubbing (rubber 'o'-rings, forceps, scissors, etc.)	When needed. Before the thermal sterilisation cycle, if required.	On clean and dry surfaces.			Avoid leaks and wipe off any excess. Do not spray near an open flame. Avoid inhaling and avoid eye area.

WE HAVE BEEN SPECIALISING WITH AIR TECHNOLOGY FOR 50 YEARS: SPECIALISATION HAS GIVEN EXCELLENT RESULTS.

To find out which unit suits your surgery's individual needs, head over to www.cattani.com.au 'Find Your Surgery Solution' page. You can work out the right plant solution for your surgery.

CATTANI TURBO SMART – INTELLIGENT SUCTION NOW IN ASIA

Cattani utility room equipment has been the beating heart of dental surgeries for 50 years, and now, finally, our suction systems and compressors are available in Asia.



Smart technology you can rely on

Based in Italy, Cattani is a family business where we believe in traditional values like loyalty and respect, so your satisfaction is very important to us. Our machinery is renowned for its performance and efficiency, and we have an outstanding reputation for service.



High quality products and service

The smooth operation of every dental practice depends on its utility room equipment, so Cattani takes the greatest care when it comes to design and manufacture, and to show how much we stand by our products, we back them with a generous warranty.



The plant room equipment of choice – for Asia

We are honored to be in Asia, and excited to bring you the air technology equipment that countless dentists around the world prefer.

HOW IS IT WE LEAD IN OUR FIELD, WHEN WE COST LESS THAN THE ALTERNATIVES? THIS IS HOW:

Constant research: this enables us to apply the latest technology to all of our products and solutions.

We enhance performance: electronic and information technology enable us to enhance the performance and reliability of our products.

We reduce costs: less maintenance and lower energy costs mean that we are always the most economical on a cost-benefit analysis.

We reduce environmental impact: we save 50% on raw materials, so that you can save between 30% and 50% on electrical consumption.



MAGNOLIA DISINFECTION FOR DENTAL SURGERIES

