

PULI-JET PLUS NEW

INSTRUCTION FOR USE

FOR SANITISING AND DISINFECTING

Puli-Jet Plus New is a non-foaming concentrate for disinfecting and simultaneous deodorization, cleaning and maintenance of dental aspirators and amalgam separators.

Puli-Jet Plus New does not contain aldehydes and phenol.

Disinfection of your suction system is vital to keep your suction lines clean and maintain your equipment and suction unit. Puli-Jet Plus New's strength as a disinfectant means that it kills bacteria, helping to maintain a safer work environment.

Puli-Jet Plus New is suitable for the following dental suction systems:

- Wet-line systems
- Semi-wet systems (Turbo Jet, Micro Smart, Turbo Smart and Maxi Smart)
- Dry-line systems (Uni Jet, Tecno Jet and Flux Jet)

Puli-Jet Plus New has a simple and accurate dispensing chamber, making it easier to measure the solution and to pour into the Pulse Cleaner Bucket. For simplicity and accuracy, Cattani recommends the use of their 'Pulse Cleaner Bucket' to automatically aspirate air and fluids to create turbulence for effective disinfection.



DISINFECTING

We recommend disinfecting the aspiration system once a week with a recommended concentrate of Puli-Jet Plus New.

Dosage: 20mL of Puli-Jet Plus New with 2.5L of water

PULI-JET PLUS NEW

INSTRUCTION FOR USE

1. Fill the Pulse Cleaner Bucket with 2.5L of warm water (approx temp 50°C).



2. Each line on the inside of the Pulse Cleaner Bucket represents 1L, so fill it up to halfway between the 2nd and 3rd line.



3. Place the Puli-Jet on a level and steady surface and slightly unscrew the red cap, but do not completely unscrew it.



4. Squeeze the bottle softly to force the solution into the measuring container at the top of the bottle, until it reaches the desired mark (mL).



5. If too much solution has entered into the container it can be poured back into the bottle by tipping the bottle towards the small tube on the side of the bottle.



6. Unscrew and remove the red cap fully and pour the solution into the Pulse Cleaner Bucket making sure that no extra solution travels up the pipe from the main bottle.



7. When all the solution has been poured into the Pulse Cleaner Bucket, replace the red cap and then with a disposable cloth, wipe away any residue from the bottle.



8. Connect your suction hoses to the Pulse Cleaner Bucket (adaptors may be required), pulsed aspiration of the solution will happen automatically.



9. Once pulsing has stopped please remove the suction hoses from the Pulse Cleaner Bucket.

Please note: there may be a small amount of solution left at the bottom of the Pulse Cleaner Bucket, this can be poured down the spitoon or aspirated manually.

It is recommended that you continue to aspirate air for several minutes before turning off the suction system. This ensures that the disinfectant reaches the suction unit in the plant room.

HOW IS IT WE LEAD IN OUR FIELD, WHEN WE COST LESS THAN THE ALTERNATIVES? THIS IS HOW:

Constant research: this enables us to apply the latest technology to all of our products and solutions.

We enhance performance: electronic and information technology enable us to enhance the performance and reliability of our products.

We reduce costs: less maintenance and lower energy costs mean that we are always the most economical on a cost-benefit analysis.

We reduce environmental impact: we save 50% on raw materials, so that you can save between 30% and 50% on electrical consumption.