

WE LOVE WHAT WE DO.

DRY SUCTION SYSTEMS

CATALOGUE




INDUSTRIAL QUALITY INDUCTION MOTORS, RENOWNED FOR THEIR RELIABILITY




WHY CATTANI

- Cattani dry separators are designed for installation close to the dental unit
- All separators are covered by a 12 MONTH WARRANTY - refer to the Cattani Warranty Statement for details.



UNI JET 75

MODEL NUMBER	020322
NO. OF SEATS	1
MAX AIR FLOW l/min	1250
MAX. VACUUM FOR CONTINUOUS OPERATION mBar	140
POWER	1ph 230V 3.1A 0.4kW 50Hz
HEIGHT MM	251
WIDTH MM	246
DEPTH MM	232
WEIGHT KG	12
SOUND OUTPUT dB(A)	62





TECNO JET

MODEL NUMBER	020590
NO. OF SEATS	2
MAX AIR FLOW l/min	2000
MAX. VACUUM FOR CONTINUOUS OPERATION mBar	140
POWER	1ph 5.5A 0.75kW
HEIGHT MM	319
WIDTH MM	303
DEPTH MM	302
WEIGHT KG	21
SOUND OUTPUT dB(A)	68



FLUX JET

MODEL NUMBER	020701
NO. OF SEATS	4
MAX AIR FLOW l/min	3300
MAX. VACUUM FOR CONTINUOUS OPERATION mBar	150
POWER	1ph 230V 7.6A 1.1kW 50Hz
HEIGHT MM	341
WIDTH MM	334
DEPTH MM	405
WEIGHT KG	25
SOUND OUTPUT dB(A)	68.5



MEDIO JET

MODEL NUMBER	021030
NO. OF SEATS	7
MAX AIR FLOW l/min	5000
MAX. VACUUM FOR CONTINUOUS OPERATION mBar	172
POWER	3ph 230/400V 9/5.2A 2.2kW 50Hz
HEIGHT MM	415
WIDTH MM	395
DEPTH MM	332
WEIGHT KG	35
SOUND OUTPUT dB(A)	76



MAXI JET

MODEL NUMBER	021090
NO. OF SEATS	12
MAX AIR FLOW l/min	8333
MAX. VACUUM FOR CONTINUOUS OPERATION mBar	166
POWER	3ph 230/400V 16.7/9.7A 4kW 50Hz
HEIGHT MM	443
WIDTH MM	566
DEPTH MM	422
WEIGHT KG	57
SOUND OUTPUT dB(A)	76.5

WE HAVE BEEN SPECIALISING WITH AIR TECHNOLOGY FOR 50 YEARS: SPECIALISATION HAS GIVEN EXCELLENT RESULTS.

To find out which unit suits your surgery's individual needs, head over to www.cattani.com.au 'Find Your Surgery Solution' page. You can work out the right plant solution for your surgery.

CATTANI TURBO SMART – INTELLIGENT SUCTION NOW IN ASIA

Cattani utility room equipment has been the beating heart of dental surgeries for 50 years, and now, finally, our suction systems and compressors are available in Asia.



Smart technology you can rely on

Based in Italy, Cattani is a family business where we believe in traditional values like loyalty and respect, so your satisfaction is very important to us. Our machinery is renowned for its performance and efficiency, and we have an outstanding reputation for service.



High quality products and service

The smooth operation of every dental practice depends on its utility room equipment, so Cattani takes the greatest care when it comes to design and manufacture, and to show how much we stand by our products, we back them with a generous warranty.



The plant room equipment of choice – for Asia

We are honored to be in Asia, and excited to bring you the air technology equipment that countless dentists around the world prefer.

HOW IS IT WE LEAD IN OUR FIELD, WHEN WE COST LESS THAN THE ALTERNATIVES? THIS IS HOW:

Constant research: this enables us to apply the latest technology to all of our products and solutions.

We enhance performance: electronic and information technology enable us to enhance the performance and reliability of our products.

We reduce costs: less maintenance and lower energy costs mean that we are always the most economical on a cost-benefit analysis.

We reduce environmental impact: we save 50% on raw materials, so that you can save between 30% and 50% on electrical consumption.