



TURBO JET **SINGLE SURGERY**

TURBO JET 1 AND COMPACT SYSTEMS



TURBO JET IS A COST EFFECTIVE SINGLE SURGERY PUMP AVAILABLE IN STANDARD (TURBO JET 1 & 2) OR COMPACT FORMAT



TURBO JET 1

This unit has a powerful motor designed for use with tip supports which are unit or cabinet mounted, compliant with ISO 10637 for vacuum.

Turbo Jet 1 motors are thermally protected and come complete with ISO4 Amalgam separation, compliant with ISO 11143.

Turbo Jet 1 can be purchased with a fitted sound-reducing cabinet, which reduces noise levels to 52 dB(A).

Please note: Sound-reducing cabinet cannot be retro fitted.



TURBO JET 1

NO. OF SEATS	1
MAX AIR FLOW l/min	650
MAX. VACUUM FOR CONTINUOUS OPERATION mBar	160
POWER	4A 0.42kW
HEIGHT MM	400
WIDTH MM	450
DEPTH MM	440
WEIGHT KG	33
SOUND OUTPUT dB(A)	61

TURBO JET 1 (NO AMALGAM SEPARATION)

NO. OF SEATS	1
MAX AIR FLOW l/min	650
MAX. VACUUM FOR CONTINUOUS OPERATION mBar	160
POWER	4A 0.42kW
HEIGHT MM	400
WIDTH MM	450
DEPTH MM	440
WEIGHT KG	33
SOUND OUTPUT dB(A)	61



TURBO JET 1 MODULAR WITHOUT BOX

NO. OF SEATS	1
MAX AIR FLOW l/min	650
MAX. VACUUM FOR CONTINUOUS OPERATION mBar	127
POWER	4A 0.42kW
HEIGHT MM	360
WIDTH MM	540
DEPTH MM	350
WEIGHT KG	25
SOUND OUTPUT dB(A)	61



TURBO JET 2 MODULAR WITHOUT BOX

NO. OF SEATS	2
MAX AIR FLOW l/min	650
MAX. VACUUM FOR CONTINUOUS OPERATION mBar	255
POWER	4A 0.42kW
HEIGHT MM	380
WIDTH MM	560
DEPTH MM	440
WEIGHT KG	30
SOUND OUTPUT dB(A)	65



TURBO JET COMPACT

Although small and light, Turbo Jet Compact is easy to use and delivers a powerful performance.

Supplied with built-in amalgam separation and prewired for installation, this suction system easily fits into a standard surgery cabinet.



TURBO JET COMPACT

NO. OF SEATS	1
MAX AIR FLOW l/min	650
MAX. VACUUM FOR CONTINUOUS OPERATION mBar	160
POWER	4A 0.42kW
HEIGHT MM	495
WIDTH MM	250
DEPTH MM	390
WEIGHT KG	33
SOUND OUTPUT dB(A)	61

WE HAVE BEEN SPECIALISING WITH AIR TECHNOLOGY FOR 50 YEARS: SPECIALISATION HAS GIVEN EXCELLENT RESULTS.

To find out which unit suits your surgery's individual needs, head over to www.cattani.com.au 'Find Your Surgery Solution' page. You can work out the right plant solution for your surgery.

CATTANI TURBO SMART – INTELLIGENT SUCTION NOW IN ASIA

Cattani utility room equipment has been the beating heart of dental surgeries for 50 years, and now, finally, our suction systems and compressors are available in Asia.



Smart technology you can rely on

Based in Italy, Cattani is a family business where we believe in traditional values like loyalty and respect, so your satisfaction is very important to us. Our machinery is renowned for its performance and efficiency, and we have an outstanding reputation for service.



High quality products and service

The smooth operation of every dental practice depends on its utility room equipment, so Cattani takes the greatest care when it comes to design and manufacture, and to show how much we stand by our products, we back them with a generous warranty.



The plant room equipment of choice – for Asia

We are honored to be in Asia, and excited to bring you the air technology equipment that countless dentists around the world prefer.

HOW IS IT WE LEAD IN OUR FIELD, WHEN WE COST LESS THAN THE ALTERNATIVES? THIS IS HOW:

Constant research: this enables us to apply the latest technology to all of our products and solutions.

We enhance performance: electronic and information technology enable us to enhance the performance and reliability of our products.

We reduce costs: less maintenance and lower energy costs mean that we are always the most economical on a cost-benefit analysis.

We reduce environmental impact: we save 50% on raw materials, so that you can save between 30% and 50% on electrical consumption.