



BLOK-JET
SILENT 2
PLANT ROOM SOLUTIONS



OUR INSTALLATION SOLUTIONS COME IN A VARIETY OF CONFIGURATIONS THAT SAVE SPACE WITHIN TIGHT CONSTRAINTS



BLOK-JET SILENT 2 INDOOR PLANT ROOM SOLUTION



The one enclosure cabinet houses suction and compressed air system for a 3-4 surgery application.

It is designed to be installed inside the clinic, with its small footprint and low noise level which makes incorporating it into a clinic as simple as possible.

The unit is fitted with lockable castor wheels for easy positioning and servicing and lockable doors for security.

External display panel shows compressor and suction status.

Simple to design into a clinic with its small footprint and low noise.

Must be installed inside the Clinic, with an ambient temperature no greater than 28oC .

Requires 2 x 10 Amp Single Phases power supplies.

Please refer to compressor and suction pages for each unit's specification data.

We aim to ensure that installation of our equipment is as simple and efficient as possible, offering a variety of solutions for even the most compact of plant rooms. Our installation solutions come in a variety of configurations that save space within tight constraints.

BLOK JET SILENT 2

NO. OF SEATS	4
HEIGHT MM	1935
WIDTH MM	760
DEPTH MM	650
WEIGHT KG	270
SOUND OUTPUT dB(A)	540



UNIT INCLUDES

AC300 & Turbo SMART B with ISO 18 Hydrocyclone



AC300

NO. OF SEATS	3-4
POWER	1ph 230V 50Hz 10.2A 1.5kW 1ph 220V 60Hz 12.2A 1.75kW
HEIGHT MM	750
WIDTH MM	620
DEPTH MM	520
WEIGHT KG	64
SOUND OUTPUT dB(A)	73.6
MAX AIR FLOW AT 5 BAR NL/min	238
TANK SIZE (LITRES)	45

* 3ph available
* w/out dryer available



TURBO SMART B

NO. OF SEATS	4
MAX AIR FLOW l/min	1750
MAX. VACUUM FOR CONTINUOUS OPERATION mBar	210
POWER	1ph 9A 1.5kW
HEIGHT MM	630
WIDTH MM	390
DEPTH MM	350
WEIGHT KG	39.5
SOUND OUTPUT dB(A)	73.7
ISO 18 HYDROCYCLONE AMALGAM RETENTION	Optional

WE HAVE BEEN SPECIALISING WITH AIR TECHNOLOGY FOR 50 YEARS: SPECIALISATION HAS GIVEN EXCELLENT RESULTS.

To find out which unit suits your surgery's individual needs, head over to www.cattani.com.au 'Find Your Surgery Solution' page. You can work out the right plant solution for your surgery.

CATTANI TURBO SMART – INTELLIGENT SUCTION NOW IN ASIA

Cattani utility room equipment has been the beating heart of dental surgeries for 50 years, and now, finally, our suction systems and compressors are available in Asia.



Smart technology you can rely on

Based in Italy, Cattani is a family business where we believe in traditional values like loyalty and respect, so your satisfaction is very important to us. Our machinery is renowned for its performance and efficiency, and we have an outstanding reputation for service.



High quality products and service

The smooth operation of every dental practice depends on its utility room equipment, so Cattani takes the greatest care when it comes to design and manufacture, and to show how much we stand by our products, we back them with a generous warranty.



The plant room equipment of choice – for Asia

We are honored to be in Asia, and excited to bring you the air technology equipment that countless dentists around the world prefer.

HOW IS IT WE LEAD IN OUR FIELD, WHEN WE COST LESS THAN THE ALTERNATIVES? THIS IS HOW:

Constant research: this enables us to apply the latest technology to all of our products and solutions.

We enhance performance: electronic and information technology enable us to enhance the performance and reliability of our products.

We reduce costs: less maintenance and lower energy costs mean that we are always the most economical on a cost-benefit analysis.

We reduce environmental impact: we save 50% on raw materials, so that you can save between 30% and 50% on electrical consumption.