

CATTANI ASPI-JET ASPIRATORS

REDUCE INFECTION RISK

If you need a simple, safe and mobile chairside aspirator, Aspi-Jet is for you.

It provides powerful, suction and air/fluid separation, and efficiently removes spray mist before it leaves the oral cavity.

INFECTION CONTROL

Maintains 300 L/min spray aspiration with all tips active.

ALL-IN-ONE SUCTION SYSTEM

Everything you need at your fingertips in a quiet, compact unit.

HIGH EFFICIENCY HEPA FILTER (OPTIONAL ADDED PROTECTION)

Reduce risk of infection with an effective antibacterial filter.



Three tip supports

Steel case powder coated white

Manual or Automatic drainage

Free Spittoon

Acoustic foam for quietness

Robust suction motor



Removes particles down to 0.01 microns and all odours with a 99.999% efficiency. (replace every 6 months)

PRODUCT SPECS

Spray aspiration:	300 L/Min
Max vacuum:	650 L/min
Vacuum (mbar):	130
Noise Level:	60 dB
Size:	H 76cm x W 50cm x D 35cm

LOW MAINTENANCE

- Daily disinfectant
- Clean internal filters
- Annual maintenance

EASY INSTALLATION

Aspi-Jet is ready to go when you are, just plug it in and turn it on.